



MICROWAVE ASSISTED EXTRACTION OF PESTICIDES FROM ENVIRONMENTAL SAMPLES



Pesticides are extensively used in modern agriculture, which could lead to serious consequences due to their biomagnification and persistent nature. Several governments require their analysis in environmental samples. Microwave assisted solvent extraction is a well-established sample preparation technique applied in several official methods. Milestone's ETHOS X equipped with fastEX-24 eT rotor was used in this study to prove its efficacy in the extraction of pesticides from environmental matrices

INTRODUCTION

Pesticides are pivotal chemical compounds for the modern agriculture production, used to control pests' diffusion. Depending on their target function, pesticides are used to treat insects, rodents, fungi and unwanted plants (weeds). For example, since 1940s, organochlorine pesticides and organophosphate pesticides are used extensively not only in agriculture but also for mosquito control. The action mechanism of these chemicals is mostly designed to disturb the physiological activities of the target organism, leading to dysfunction and reduced vitality. Despite their fundamental use in the modern agriculture, these molecules can cause neurological damage, endocrine disorders, and have acute and chronic health effects on human¹. Moreover, some of these chemicals belong to the class of persistent

organic pollutants (POPs) with high persistence in the environment².

EPA 3546³ outlines the procedure for extracting water insoluble or slightly water-soluble organic compounds from soils, clays, sediments, sludges, and solid wastes.

EPA 3546 is a specific method for Microwave Assisted Solvent Extraction (MASE), a well-established sample preparation technique that enables extractions with reduced solvent volume and time. This application note represents a guideline for the extraction of the priority pesticides from both standard reference materials and spiked materials using the official method EPA 3546.

APPLICATION REPORT

ETHOS X – Pesticides – EPA 3546 Method



EXPERIMENTAL

EQUIPMENT

- Milestone's ETHOS X.
- fastEX-24 eT rotor.⁴
- 100-mL disposable glass vials.
- SFS-24 (Simultaneous Filtration System).
- GC-MS/MS
- HPLC MS/MS




Figure 1 – Milestone ETHOS X with fastex-24 eT (left) and SFS-24 filtration system (right).

STANDARD AND REAGENTS

Standards, surrogates and internal standard were purchased by Sigma Aldrich. Grade solvent pesticide were used. Sodium sulfate anhydrous, silica gel (activated for at least 16 h at 130°C) and glass wool or paper filter were used in the clean-up procedure. According to the analytical method EPA 8270e⁵, internal surrogates and standards were used.

Analyte	CAS-No
α -BHC	319-84-6
γ -BHC	58-89-9
β -BHC	319-85-7
δ -BHC	319-86-8
Aldrin	309-00-2
Heptachlor	76-44-8
γ - α -chlordane	5103-74-2 \ 57-74-9

α -Endosulfan	1031-07-8
4,4'-DDE	72-55-9
Dieldrin	60-57-1
Endrin	72-20-8
β - Endosulfan II	33213-65-9
4,4'-DDD	72-54-8
2,4'-DDT	789-02-6
4,4'-DDT	50-29-3
Mirex	2385-85-5
Dichlorvos	62-73-7
Endosulfan Sulfate	1031-07-8
Endrin aldeide	7421-93-4
Endrin ketone	7378-10-1
Heptachlor Epoxide	1024-57-3
Methoxychlor	72-43-5
Demeton	8065-48-3
Dimethoate	60-51-5
Malathion	121-75-5
Parathion	56-38-2
Parathion-methyl	298-00-0
Mevinphos	7786-34-7
Phorate	298-02-2
Fenitrothion	122-14-5
Isocarbophos	24353-61-5
Methidathion	950-37-8
Toxaphene	8001-35-2

Table 1 - Pesticides Stock solution

APPLICATION REPORT

ETHOS X – Pesticides – EPA 3546 Method



Analyte	CAS-No
Isoproturon-d6	217487-17-7
Biphenyl-d10	1486-01-7
Atrazina-d5	163165-75-1
Phenanthrene-d10	1517-22-2
Pirimicarb-d6	1015854-66-6
PCB 138	35065-28-2
Triphenylphosphate	115-86-6

Table 2- Internal Standard Solution

SAMPLES

The clay loam 1 - CRM847 certified reference material were used for the determination of pesticides. For pesticides not included in the certified materials a spiking stock solution on blank soil was used.

SAMPLE PREPARATION

The samples were collected and stored in accordance with the requirements of EPA 3546. Decant and discard any water layer on a sediment sample. Discard any foreign objects such as sticks, leaves, and rocks. Mix the sample thoroughly, especially composited samples. Grind or otherwise reduce the particle size of the waste so that it either passes through a 1 mm sieve or can be extruded through a 1 mm hole.

Ground samples, wet or dried, were weighed directly into the 100 mL extraction disposable glass vials of the fastEX-24 eT rotor. 30 ml of acetone-hexane (1:1) was used as extraction mixture. An aliquot of the internal standard solution was added to the samples just prior to solvent addition then the glass vials were closed (automatic capping tool available).

EXTRACTION PROCESS AND CLEAN UP

According to the moisture content, the proper built-in method was selected.

Step	Time (min)	Power (W)	Temperature (°C)
1	15	up to 1600	110
2	10	up to 1600	110

Table 3 - Microwave Program




Figure 2 - Microwave run profile

After the extraction, samples were filtered with milestone SFS-24 simultaneous filtration system using sodium sulfate anhydrous. The vials were rinsed with additional solvent aliquots. SFS-24 allows to filter 24 samples simultaneously with different types of filters available. Extracts and rinse solution were collected together.

The extract was subsequently concentrated with nitrogen flow. If purification is not required, concentrate directly until 0.5 ml and add the appropriate surrogate standard solution to achieve the surrogate standard concentration. If purification is necessary, concentrate the extract directly until 2 ml. Purify the solution according to the method (EPA 3610, 3620, 3630, 3640, 3660). Finally, the extracts obtained by ETHOS X were concentrated for analysis.

APPLICATION REPORT

ETHOS X – Pesticides – EPA 3546 Method



ANALITICAL CONDITIONS

Based on the pesticides compositions both GC-MS/MS and HPLC MS/MS with triple quadrupole were used.

GC-MS/MS is equipped with a split-splitless injector, autosampler and triple quadrupole mass spectrometer. Sample injection volume was 1 µl. A 30 m x 0.25 x 0.25 RXI 5MS Capillary column (Restek) was used for the analyses. A five steps ramp oven program was used:

Rate (°C/min)	Temperature (°C)	Plateaus (min)
20	50	2
30	150	0
7	260	0
20	290	11

Table 4 – GC-MS/MS oven program

Helium was used as the carrier gas at a linearity velocity of approximately 65 cm/s.

HPLC MS/MS was used to ensure the recovery data and to better quantify some compounds. An oven temperature of 40 °C is used with an injection volume of 0.25 µl.

Mobile phase:

- Water 0,1% formic acid 5 mM ammonium formate.
- Methanol 0,1 % formic acid 5 mM ammonium formate.

Time (min)	Water (%)	Methanol (%)
0.2	95	5
11	0	100
13	0	100
13.1	95	5

Table 5 – HPLC-MS/MS gradient method

RESULTS AND DISCUSSION

Results from extractions of Clay loam 1 - CRM847 are shown in Table 6. Recovery for all compounds are in the range 70-120% of the certified standard reference material.

Analyte	Certified value (µg/kg)	Ethos X (µg/kg)	Recovery (%)	RSD (%)
δ-BHC	138	128.89	93.4	8.7
α-BHC	221	188.95	85.5	7.0
β-BHC	295	342.2	116	2.2
α-Chlordane	309	330.63	107	3.3
γ-Chlordane	171	182.28	106.6	2.4
4,4'-DDD	120	107.28	89.4	3.1
4,4'-DDE	315	290.74	92.3	3.3
4,4'-DDT	92.1	84.27	91.5	11.5
Dieldrin	53.3	54.79	102.8	9.6
Endosulfan I	211	162.25	76.9	4.8
Endosulfan II	225	160.65	71.4	7.9
Endosulfan Sulfate	159	159.63	100.4	4.0
Endrin Ketone	170	139.57	82.1	6.3
Endrin	162	156.81	96.8	4.4
Heptachlor epoxide	127	119.25	93.9	3.7
Methoxychlor	290	260.13	89.7	6.6

Table 6 - Pesticides recovery from Clay loam 1 - CRM847 (1g) (n=4).

APPLICATION REPORT

ETHOS X – Pesticides – EPA 3546 Method



Additionally, a pesticide mixture was spiked to a blank soil in order to test the performance of the fastEX-24 eT on a wider list of pesticides (Table 7).

Analyte	Spike concentration (µg/kg)	Ethos X (µg/kg)	Recovery (%)	RSD (%)
α-BHC	50	52.00	104.0	4.6
γ-BHC	50	50.00	100.0	8.1
β-BHC	50	60.00	120.0	3.2
δ-BHC	50	56.67	113.3	4.1
Aldrin	50	53.33	106.7	5.2
Heptachlor	50	56.67	113.3	2.6
alfa clordano	50	48.40	96.8	4.8
gamma clordano	50	52.67	105.3	6.1
4,4'-DDE	50	56.67	113.3	3.2
Dieldrin	50	53.33	106.7	1.3
Endrin	50	52.05	104.1	1.9
β- Endosulfan II	50	50.00	100.0	4.1
4,4'-DDD	50	56.67	113.3	3.8
2,4'-DDT	50	46.65	93.3	2.4
4,4'-DDT	50	45.70	91.4	2.9
Mirex	50	53.05	106.1	5.1
Dichlorvos	50	40.00	80.0	9.3
Demeton	50	41.45	82.9	6.8
Dimethoate	50	52.50	105.0	4.2
Malathion	50	45.70	91.4	3.6
Parathion	50	43.95	87.9	2.7
Parathion-methyl	50	53.33	106.7	2.9
Mevinphos	50	54.20	108.4	3.9
Phorate	50	46.67	93.3	4.6
Fenitrothion	50	56.67	113.3	5.8
Isocarbophos	50	46.25	92.5	6.4
Methidathion	50	43.15	86.3	8.9
Endosulfan Sulfate	50	46.90	93.8	10.6
Endrin aldeide	50	39.45	78.9	3.5
Endrin ketone	50	39.70	79.4	4.9
Heptachlor Epoxide	50	47.30	94.6	6.3
Methoxychlor	50	53.35	106.7	5.5
Toxaphene	50	41.80	83.6	4.9

Table 7 - Pesticides recovery from Spike solution (n=4)

APPLICATION REPORT

ETHOS X – Pesticides – EPA 3546 Method



CONCLUSION

The results demonstrate the efficiency of the ETHOS X with fastEX-24 eT rotor for the pesticides extraction from environmental matrices. High recovery rate for all the tested molecules showed the great extraction efficiency.

The fastEX-24 eT enables simultaneous solvent extraction of up to 24 samples in only 40 minutes (cooling step included). In turns this means that is able to extract over 200 samples in 8-hour workday. Contamination, memory effects, and cleaning are completely eliminated due to the use of disposable glass vials. The use of contactless temperature control ensures high reproducibility and full recovery of the target analytes for full compliance with Official Methods.

Thanks to the unique design, fastEX-24 eT is easily applied to even more challenging matrices such as solid wastes and plastics. ETHOS X provides extracts with the lowest solvent usage and significant time saving compared to all the other extraction techniques.

The ETHOS X with all its unique features fully addresses the need of environmental laboratories in terms of productivity, ease of use, running costs, and extraction quality.

REFERENCES

- 1- Chemical Pesticides and Human Health: The Urgent Need for a New Concept in Agriculture
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4947579/>
- 2- Stockholm Convention – POPs project
<http://www.pops.int/>
- 3- EPA 3546- Microwave extraction
<https://www.epa.gov/sites/production/files/2015-12/documents/3546.pdf>
- 4- ETHOS X and fastEX 24 eT
<https://www.milestonesrl.com/products/microwave-extraction/ethos-x-for-environmental>
- 5- EPA 8270 E – Semivolatile organic compound GC-MS
<https://www.epa.gov/esam/epa-method-8270e-sw-846-semivolatile-organic-compounds-gas-chromatographymass-spectrometry-gc>

ABOUT MILESTONE

At Milestone we help chemists by providing the most innovative technology for metals analysis, direct mercury analysis and the application of microwave technology to extraction, ashing and synthesis. Since 1988 Milestone has helped chemists in their work to enhance food, pharmaceutical and consumer product safety, and to improve our world by controlling pollutants in the environment.



MILESTONE
H E L P I N G
C H E M I S T S

MILESTONE Srl - Via Fatebenefratelli, 1/5 - 24010 Sorisole (BG) - Italy
Tel: +39 035 573857 - Fax: +39 035 575498
www.milestonesrl.com email: analytical@milestonesrl.com