

## OPTIMIZED ACID DIGESTION OF SILICATE ROCKS FOR HIGH-PRECISION GEOCHEMICAL ANALYSIS

**Microwave-assisted acid digestion of geological samples for elemental quantification on ICP-MS: silicate rocks without resistant minerals.**

### | SUMMARY

Historically, the preferred analytical technique for mineral analysis was arc emission spectroscopy, but this is no longer suitable in the isotope ratios determination. This is one of the reasons why, in contemporary laboratories, more advanced techniques such as Inductively Coupled Plasma Optical Emission Spectroscopy (ICP-OES) and Inductively Coupled Plasma Mass Spectrometry (ICP-MS) have become standard. These techniques, particularly for trace element analysis, necessitate complete sample dissolution—a critical step that is typically achieved through strong acid digestion.

It is important to make an analytical distinction when dealing with samples containing refractory minerals that are not

readily dissolved in open-vessel HF-HNO<sub>3</sub> mixtures (e.g., ZrSiO<sub>4</sub>, FeCr<sub>2</sub>O<sub>4</sub>, BaSO<sub>4</sub>). For such materials, closed-vessel digestion is required to reach significantly higher temperatures and pressures and effectively digest the sample. This is where Milestone's Single Reaction Chamber (SRC) technology becomes essential. The SRC system utilizes microwave-assisted digestion in a fully stainless-steel, 1-liter pressurized reactor. With the UltraWAVE 3 platform, temperatures up to 300 °C and pressures up to 199 bar can be achieved.

In collaboration with Petrology Professor Kamber Balz (*Queensland University of Technology*) and with contributions from geological laboratories, Milestone has developed a comprehensive sample



preparation method applicable to a wide range of geological materials. These have been classified into five categories based on their chemical behavior and composition. This document presents a four-step protocol specifically designed for the complete dissolution of silicate rocks samples that do not contain refractory minerals (Basalts, Greenschists, Greywackes), ensuring complete digestion of up to 20 samples within 24 hours, thus significantly reducing the time required for elemental analysis sample preparation.

## I EXPERIMENTAL

### INSTRUMENTATION

- UW3 system
- 20-position rack
- PTFE 15 mL vials with caps
- Chiller 1 kW
- Nitrogen gas line (40 bar / or min. 10 bar using gas booster)
- Hotblock with fumehood connection
- ICP-MS

### REAGENTS

- HF 48%
- HNO<sub>3</sub> 67%
- HNO<sub>3</sub> 50%
- Distilled water

## I REAGENT HANDLING

Reagent addition, a common task during the sample preparation process, poses safety concerns for the operator, is time-consuming, and can lead to contamination if not performed correctly. The use of an

automatic dosing station, such as Milestone easyFILL, mitigates and limits these risks, especially when the operator is exposed to toxic acid as HF, very common on geochemical methods. Specifically designed to precisely add reagents to digestion vessels and vials, easyFILL helps chemists optimize their procedures by reducing the risk of human error in trace analysis. Capable of dispensing various types of reagents into digestion vessels and vials, easyFILL minimizes manual handling of the digestion mixture, thereby reducing the risk of contamination. EasyFILL is fully compatible with all UW3 racks, including the 20-position rack.



• Figure 2: Milestone's easyFILL

## SINGLE REACTION CHAMBER (SRC) TECHNOLOGY: UW3 OVERVIEW



Figure 1 UW3 microwave digestion system with focus on the SRC chamber



SRC technology represents the latest revolution in microwave closed-vessel digestion. Unlike traditional microwave ovens, SRC utilizes a pressurized one-liter stainless steel microwave reactor which is also the digestion vessel. The process is straightforward: samples and reagents are loaded into vials and placed inside the reactor cavity together with a suitable baseload that allows heat homogeneity around the samples. Once loaded, the reactor is securely closed and sealed. Automatic introduction of nitrogen gas (40 bar) follows for internal pressurization, which also serves to raise the boiling points of the solutions and not lose the volatile elements. At the end of the process, the reactor is cooled using a powerful water-cooling system.

#### METHOD DESCRIPTION

The acid digestion protocol developed for the preparation of geological samples, specifically, siliceous rocks not containing refractory minerals, comprises four essential and sequential steps, each designed to ensure complete dissolution and compatibility with subsequent elemental or isotopic analysis and avoiding the use of  $\text{HClO}_4$ .

The **first stage** involves a high-temperature digestion aimed at the decomposition of silicate structures. This is achieved through the addition of concentrated hydrofluoric acid to the samples placed within a 20-position rack. Digestion is performed using the SRC technology, which maintains a stable temperature of 250 °C for an

extended duration. Under these conditions, silicon is effectively complexed with fluoride ions, as hydrofluoric acid is capable of cleaving strong metal–Si–O bonds.

The **second step** consists of an evaporation phase carried out on a hotblock, in which the samples are brought to dryness; This step does not affect the recovery of elements, not even volatile ones. This step promotes only the volatilization and removal of silicon in the form of silicon tetrafluoride ( $\text{SiF}_4$ ), effectively reducing the sample matrix and facilitating the digestion of the remaining elements.

The **third step** involves the addition of 2 mL of nitric acid to each sample, followed by gentle heating on the hotblock. This step serves to reconvert any residual fluoride species into soluble complexes and ensures the breakdown of intermediate fluoride salts, thereby preparing the sample for the final digestion.

The fourth and **final step** consists of a second digestion performed in the UltraWAVE 3 system. This advanced SRC-based microwave digestion unit is capable of fully dissolve compounds such as aluminum fluoride ( $\text{AlF}_3$ ) and other persistent residues that may have formed during earlier steps. The application of this final digestion step ensures complete mineral decomposition, producing a homogeneous solution suitable for trace



elemental determination via ICP-OES or ICP-MS.




Figure 3: Schematic workflow for silicate rock digestion and conversion to weak HNO<sub>3</sub> solution.

### OPERATING CONDITIONS ON UW3

For the first digestion, a mixture of 3 mL of hydrofluoric acid and 1.5 mL of nitric acid is added to 0.100 g of sample in each 15mL test tube. The digestion program followed by the instrument is as follows:

First digestion step, operating conditions:

Step	Time	Power (W)	Temp T1 (°C)	Temp T2 (°C)	Pressure (bar)
1	00:35:00	1500	250	70	120
2	00:60:00	1500	250	70	140

The second digestion step, performed using SRC technology and representing the final stage of the overall process, involves the addition of 3 mL of HNO<sub>3</sub> and 6 mL of deionized water, followed by the following program:

Last digestion step, operating conditions:

Step	Time	Power (W)	Temp T1 (°C)	Temp T2 (°C)	Pressure (bar)
1	00:35:00	1500	200	70	120
2	00:60:00	1500	200	70	120

## RESULTS AND DISCUSSION

Table 1: Recommended ReoREM values (Jochum et al., 2016) except: Zr from Kamber et al. (2005), Ni, Cd, Sb, Pb from Marx and Kamber (2010).

Element	Reference value (µg/g)	Uncertainty (µg/g)	Silicate method digestion (n=4) (µg/g)	RSD (µg/g)	Accuracy (fraction)
Li	9.13	0.22	9.28	0.16	1.016
Be	2.17	0.10	2.23	0.06	1.028
Ti	13575.08	143.84	13475.64	167.20	0.993
Cr	15.85	0.38	15.81	0.08	0.997
Ni	11.96	0.24	11.92	0.16	0.997
Cu	19.66	0.72	21.22	0.23	1.079
As	0.860	0.20	0.967	0.022	1.124



Sr	337.4	6.70	339.66	5.35	1.007
Zr	185	1.68	179.82	2.62	0.972
Mo	250.6	6.70	252.5	3.07	1.008
Cd	0.227	0.022	0.211	0.01	0.931
Sb	0.274	0.03	0.27	0.01	0.983
Ba	683.9	4.70	688.9	9.50	1.007
Ce	53.12	0.33	53.58	0.77	1.009
Eu	1.989	0.024	1.962	0.031	0.987
Pb	9.732	0.272	9.852	0.140	1.012
U	1.683	0.017	1.689	0.023	1.004

Table 1 lists the main elements analyzed by ICP-MS following acid digestion according to the method presented in this document. For a complete list of analyzed elements, please refer to the eBook (*Efficient rock digest preparation for geochemists: a practical handbook*). All recoveries were accurate, exceeding 90%, with an average of 99%.




Figure 4: Graph BCR-2, elements sorted by atomic number and relative recovery.

Figure 4 graphically displays the recoveries of each element, ordered by atomic number. Notably, even for volatile elements (Cd, Pb), recovery rates fall within the optimal range, yielding reproducible results.



### I CONCLUSION

SRC technology successfully achieved complete digestion of siliceous rock samples, obtaining excellent recoveries of all analyzed elements and optimizing sample preparation times without losing any performance quality.

Using Milestone easyFILL system, reagent handling has been improved, significantly reducing the risk of operator exposure to hazardous acids (such as HF) and automating their addition in all 20 positions.

### REFERENCES

[1] *Efficient rock digest preparation for geochemists: a practical handbook*, Milestone Helping Chemists, Balz Kamber, Diego Carnaroglio, Giulio Colnaghi, Mirco Rossetti, Gianpaolo Rota, Osama Ghidan.

